English

User Manual **Weigh Modules** WMC

Deutsch

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Français

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1 Introduction

1.1 This User Manual

This User manual contains all information for the **operator** of the product.

- Read this User manual carefully before use.
- Keep this User manual for future reference.
- Pass this User manual to any further owner or user of the product.

1.2 Further documents

In addition to this printed User manual you can download the following documents from

Documentation WMC

► www.mt.com/ind-wmc-support

- Technical data sheet
- Installation information (for trained personnel under the control of the operating company)
- Reference manual for the command set

1.3 Manufacturer information

The contact information of the manufacturer of the product is as follows:

- **Name:** Mettler-Toledo GmbH
- **Weblink:** <http://www.mt.com>
- **Physical address:** Im Langacher 44, 8606 Greifensee, Switzerland

2 Safety Information

2.1 Intended use

- Use the product only for weighing in accordance with this User manual.
- The weigh module is intended for indoor use only.
- Any other type of use and operation beyond the limits of technical specifications is considered as not intended.

2.2 Definition of signal words and warning symbols

Safety notes are indicated by signal words and warning symbols and contain warnings and information about safety issues. Ignoring safety notes can lead to personal injury, damage to the instrument, malfunctions and erroneous results.

Signal words

WARNING

for a hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.

CAUTION

for a hazardous situation with low risk, resulting in minor or moderate injury if not avoided.

NOTICE

for a hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

Note

(no symbol)
for useful information about the product.

Warning symbols

General hazard

Electrical shock

2.3 Product specific safety notes

Your weigh module meets the state of the art technology and complies with all recognized safety rules, however, certain hazards can arise.

Do not open the weigh module: It does not contain any parts which can be maintained, repaired or replaced by the user. If you ever have problems with your weigh module, contact your authorized METTLER TOLEDO dealer or service representative.

Observe instructions

Always operate and use your weigh module only in accordance with the instructions contained in the product documentation. The instructions for setting up your weigh module must be strictly observed.

If the weigh module is not used according to the product manuals, protection of the weigh module may be impaired and METTLER TOLEDO assumes no liability.

Staff safety

Use only METTLER TOLEDO accessories and peripheral devices, these items are designed to work optimally with your weigh module.

Explosion hazard

It is not permitted to use the weigh module in explosive atmospheres of gases, steam, fog, dust and flammable dust (hazardous environments).

Safety notes

CAUTION

Risk of electric shock

The weigh modules may only be connected to DC power sources that meet the 12VDC \pm 3% at all times.

The power supply must be approved by the respective national test center of the country in which the weigh module will be used.

3 WMC Weigh Modules

3.1 Specifications

Parameter		WMC25-SH	WMC24-SH	WMC15-SH
Maximum capacity		21 g	21 g	11 g
Readability		0.01 mg	0.1 mg	0.01 mg
Measurement properties				
Specification temperatures		10 ... 30 °C		
Specification humidity		20 ... 80% rH		
Limit values				
Repeatability ¹⁾ (at nominal load)	sd	0.02 mg (20 g)	0.1 mg (20 g)	0.02 mg (10 g)
Repeatability (at low load)	sd	0.012 mg (5 g)	0.08 mg (5 g)	0.012 mg (5 g)
Linearity deviation		+/- 0.05 mg	± 0.2 mg	± 0.05 mg
Eccentricity deviation 2 mm offset		0.12 mg (10 g)	0.4 mg (10 g)	0.12 mg (5 g)
Sensitivity offset (test load)		$10 \times 10^{-6} \cdot R_{nt}$	$20 \times 10^{-6} \cdot R_{nt}$	$15 \times 10^{-6} \cdot R_{nt}$
Sensitivity temperature drift ²⁾		$0.0002\% / ^\circ\text{C} \cdot R_{nt}$	$0.0002\% / ^\circ\text{C} \cdot R_{nt}$	$0.0002\% / ^\circ\text{C} \cdot R_{nt}$
Sensitivity stability		$0.0001\% / \alpha \cdot R_{nt}$	$0.0001\% / \alpha \cdot R_{nt}$	$0.0001\% / \alpha \cdot R_{nt}$
Dynamics				
Settling time ³⁾ 3 × sd		0.7 s	0.6 s	0.6 s
Interface update rate (max.)		23 weight values per second (with terminal), 92 weight values per second (without terminal)		

sd = Standard deviation

R_n = Net weight (sample weight)

t

R_g = Gross weight

α = Year (annum)

r

¹⁾ Applies only after adjustment at nominal capacity with an OIML E2 weight.

²⁾ Temperature range 10 ... 30 °C.

³⁾ The settling time is the time between when the object to be weighed is placed on the scale and a stable signal is transmitted – this assumes optimal ambient conditions and optimal parameter settings.

3.2 Unpacking the weigh module

The weigh unit and the electronic unit are each packed in separate cardboard boxes along with the associated hardware and appropriate mounting materials. The boxes are labeled with their contents as follows: "Weighing module," and "Electronic Unit". Both cardboard boxes are placed inside a large transport box. The terminal must be ordered separately, and consequently it is supplied separately.

Unpacking the electronic unit and terminal:

Remove the top foam cushion, then remove the unit and the hardware.

Unpacking the weigh unit:

Remove the grey foam cushion. Carefully remove the weigh unit from the cushion and place it smoothly on a level surface.

Note

We recommend to keep the original packaging and to use it to transport the balance. The weighing pan must always be removed before transportation. Mount the protection cap to protect the interior of the weigh module.

3.3 Scope of delivery**“Weighing Module” Box**

Parts	WMC24-SH	WMC15-SH	WMC25-SH
WMC weigh unit	✓	✓	✓
Standard weighing pan	✓	✓	✓
Weighing pan adapter	✓	✓	✓
Draft shield	✓	✓	✓
WMC User Manual (this document)	✓	✓	✓
Production certificate and CE Declaration of Conformity	✓	✓	✓

“Electronic Unit” Box

Parts	WMC24-SH	WMC15-SH	WMC25-SH
Electronic unit	✓	✓	✓
Mounting bracket for electronic unit, including clip and screws for attaching to DIN standardized rail	✓	✓	✓
AC adapter	✓	✓	✓
Power cable (country-specific)	✓	✓	✓

Terminal option: “Terminal” Box

SWT terminal (monochrome display), including protective cover
2 terminal-electronic unit connector cables (0.575 m/1.9 ft and 2 m/6.5 ft long)

4 Operation

4.1 Electrical connections

WMC weigh modules have a standard RS232C interface for the communication with the control systems. Additionally, a second optional interface can be added, which can be either a second RS232C interface or an Ethernet TCP/IP interface.

Interface descriptions

Interface type	RS232C, bidirectional, full duplex, 600 to 38'400 ¹⁾ bps (baud rate can be selected using interface commands)
Interface update rate (max.)	23 weight values per second (with terminal) 92 weight values per second (without terminal)

¹⁾ 38400 baud only possible:

- Weigh module without terminal, or
- Weigh module with terminal, only via the optional RS232C interface.

Power supply for the weigh module

Input voltage	12 V DC (nominal)
Input current	2 ADC

- Use a stable power supply with no voltage fluctuations.
- If voltage fluctuations cannot be prevented, use a voltage regulator to deliver a constant voltage value to the load cell.
- The power supply must be approved by the respective national test center of the country in which the weigh module will be used.

Always operate the weigh module with the standard AC adapter delivered in the scope of delivery. The AC adapter is suitable for use **only** with the following input characteristics:

100 – 240 V AC, -15%/+10%, 50/60 Hz, 0.5A

Note

Check whether your local power supply falls within this range. If it does not, DO NOT connect the electronic unit or the AC adapter to your power supply and contact your local METTLER TOLEDO office.

If you use an alternative AC adapter other than the one in the standard scope of delivery, make sure that it is a tested AC adapter with SELV output current and ensure correct polarity.

4.2 Environmental conditions

WMC weigh modules can be operated within the following environmental conditions:

Temperature Range	Operating / Compensated	+10 ... +30 °C
	Allowable ambient	+5 ... +40 °C
Relative air humidity		Max. 80 % at 31 °C, decreasing linearly to 50 % at 40 °C, non-condensing
Height above mean sea level		Max. 4'000 m (13'330 ft)
Warm-up time		At least 60 minutes after power up At least 60 minutes after the weigh module has been connected to the power supply; the weigh module can be used immediately if it is turned on from standby mode (when operated with a terminal connected).

4.3 Ingress protection

IP rating of the WMC weigh modules:

In operation, the weigh unit is rated IP30. If the wash-down module is mounted on the top (available as an accessory), an enhanced IP rating of IP54 can be reached during operation.

The electronic unit meets IP40.

The SWT terminal is rated IP54.

If the wash-down module is mounted on the top (available as an accessory), an enhanced IP rating of IP56 can be reached in wash-down configuration (weighing not possible).

5 Maintenance

To ensure that your weighing module remains reliable, accurate, and functional for many years to come, the individual components must be cleaned and maintained periodically as appropriate for the intensity of use and the risk of contamination with debris.

5.1 Cleaning the weigh module

Clean the weighing pan and the weigh unit housing regularly with a damp cloth. The electronic unit and terminal can also be cleaned in this fashion as needed. For tougher dirt, a mild household cleaner may be used. Make sure that no liquid penetrates inside the components (protect with plastic cover)!

With the wash-down module (available as accessory), the IP rating can be increased to IP56 for cleaning. In this configuration, the weigh module can be cleaned with a water jet.

Please observe the following notes:

⚠ WARNING

Risk of electric shock

- 1 Disconnect the weigh module from the power supply prior to cleaning and maintenance.
- 2 Only use METTLER TOLEDO power cable, if these need to be replaced.
- 3 Ensure that no liquid comes into contact with the weigh module, terminal or AC adapter.
- 4 Do not open the weigh module, terminal or AC adapter.
These contain no user-serviceable parts.

⚠ CAUTION

Damage to weigh module

Under no circumstances use cleaning agents containing solvents or abrasive agents, as this can damage or scratch certain surfaces (particularly the terminal screen).

5.2 Maintenance

Your weigh module is a precision instrument, and periodic maintenance is one of the basic requirements to ensure it will perform well for you for many years to come.

Maintenance intervals will depend on the duration of use and the application and ambient conditions. Maintenance must be performed by a technician trained by METTLER TOLEDO.

Ask your METTLER TOLEDO office about service packages – regular maintenance by an authorized service technician will ensure your weigh module remains accurate and lengthen its service life.

5.3 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

Thank you for your contribution to environmental protection.

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© Mettler-Toledo GmbH 06/2017